

## OPHTHALMIC GAS



SF<sub>6</sub> Kit includes 1 50 ml pure sulphur hexafluoride cylinder for long-term replacement of the vitreous humour extracted during vitreoretinal surgery and 4 gas injection sets. Product can be used until cylinder content is consumed or as a single use product by those who favour disposable products. Accessories permit to use a fresh sterile gas injection set for each gas withdrawal. Special 50 ml cylinder with multi-layer bag resists pressure variations, ensuring shipping by air without any problems or limitations.

	GOT SF <sub>6</sub> Multi		GIS
	Gas cylinder	Injection set	Gas Injection Set
<b>Ingredient</b>	100 % sulphur hexafluoride	plastic	plastic
<b>Purity</b>	> 99%	-	-
<b>Form</b>	gas	solid	solid
<b>Amount</b>	50 ml	1 set	3 sets
<b>Packaging</b>	aluminium cylinder with multilayer bag	pouch	pouch
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>▶ sterile cylinder by ethylene oxide (according to EN ISO 11135-1:2008)</li> <li>▶ non sterile gas</li> </ul>	ethylene oxide (according to EN ISO 11135-1:2008)	ethylene oxide (according to EN ISO 11135-1:2008)
<b>Storage temperature</b>	room temperature		
<b>Expiry date</b>	24 months		
<b>Intended use</b>	intraocular tamponade during vitreoretinal surgery	set for intraocular injection of ophthalmic gas	set for intraocular injection of ophthalmic gas
<b>Tamponade time</b>	10 – 14 days		
<b>Non-expansile concentration</b>	20%		
<b>Gas volume for non-expansile concentration</b>	10 ml		
<b>Class</b>	medical device		
<b>GMDN code</b>	P 45125		P 47242
<b>MD registration number</b>	282		
<b>Medical device class</b>	IIb		IIa
<b>CE mark</b>	CE 0546		CE 0546

**Ref.:** SFP 001**Contents:**

- ▶ 1 cylinder of 50 ml pure gas
- ▶ 4 gas injection sets each made of 1 50-ml polypropylene syringe, 1 0.2 µm filter, 1 connector,
- ▶ 4 30G needles,
- ▶ 4 patient wristbands with warnings,
- ▶ 12 labels for product traceability,
- ▶ 2 inserts

**EAN 13 barcode:**

8021118000388

ALCHIMIA

Viale Austria, 14  
35020 Ponte S. Nicolò (PD)  
Italy

Tel. +39 049 8962074 / 64  
Fax +39 049 8962071

www.alchimiasrl.com  
info@alchimiasrl.com