

HIGH PURITY STERILE
SILICONE OIL, 5000 cS,
15 ML FOR INTRAOCULAR
TAMPONADE

RS-OIL 5000

SILICONE OIL



RS-OIL is highly purified silicone oil (100% pure polydimethylsiloxane), 5000 cS, 15 ml, sterile and apyrogenic, packaged in vial, intended for long-term replacement of the vitreous humour extracted during vitreoretinal surgery. Product warrants a stable tamponade over the time and reduced incidence of complications.

RS-OIL 5000	
Ingredient	100% silicone oil (100% polydimethylsiloxane)
Viscosity (25°C)	5000 cS
Specific gravity (25°C)	0,967 - 0,975
Refractive index (25°C)	1,40
Volatility (1 g/2 h/150°C)	< 0,1%
Form	fluid
Volume	15 ml
Packaging	vial in medical pouch
Sterilization	moist heat according to EN ISO 17665-1:2006
Storage temperature	room temperature
Expiry date	36 months
Intended use	intraocular tamponade during vitreoretinal surgery
Class	medical device
GMDN code	P 45125
MD registration number	3534/R
Medical device class	IIb
CE mark	CE 0546

Ref.: RSO 004-00

Contents:

- ▶ 1 RS-OIL vial (100% high purity, sterile Silicone Oil 5000 cS, 15 ml)
- ▶ 4 labels for product traceability
- ▶ 1 insert

EAN 13 barcode:

0211180000676

ALCHIMIA

Viale Austria, 14
35020 Ponte S.Nicolò (PD)
Italy

Tel. +39 049 8962074 / 64
Fax +39 049 8962071

www.alchimiasrl.com
info@alchimiasrl.com