

HIGH PURITY STERILE
SILICONE OIL 5000 cS,
10 ML IN SYRINGE
WITH ADAPTER

RS-OIL ECS 5000

EASY CONNECTION SYRINGE

SILICONE OIL



RS-OIL ECS 5000	
Ingredient	100% silicone oil (100% polydimethylsiloxane)
Viscosity (25°C)	5000 cS
Specific gravity (25°C)	0,967 - 0,975
Refractive index (25°C)	1,40
Volatility (1 g/2 h/150°C)	< 0,1%
Form	fluid
Volume	10 ml
Packaging	chemically inert unbreakable syringe with adapter in medical pouch
Sterilization	moist heat according to EN ISO 17665-1:2006
Storage temperature	room temperature
Expiry date	24 months
Intended use	intraocular tamponade during vitreoretinal surgery
Class	medical device
GMDN code	P 45125
MD registration number	3607/R
Medical device class	IIb
CE mark	CE 0546

RS-OIL ECS is highly purified silicone oil (100% pure polydimethylsiloxane), 5000 cS, 10 ml, sterile and apyrogenic, packaged in pre-filled syringe, intended for long-term replacement of the vitreous humour extracted during vitreoretinal surgery. Product warrants a stable tamponade over the time and reduced incidence of complications. The ready-to-use syringe is made of inert, plastic and silicone oil compatible material to increase patient's safety by reducing the risk of breakage during handling or surgery.

Ref.: RSO 007-00

Contents:

- ▶ 1 chemically inert unbreakable syringe of high purity sterile RS-OIL ECS (100% Silicone Oil 5000 cS), 10 ml
- ▶ 4 labels for product traceability
- ▶ 1 insert

EAN 13 barcode:

8021118000098

ALCHIMIA

Viale Austria, 14
35020 Ponte S.Nicolò (PD)
Italy

Tel. +39 049 8962074 / 64
Fax +39 049 8962071

www.alchimiasrl.com
info@alchimiasrl.com